



SIMPLIEV VOLIR WORK LISING A FILLLY ALITOMATED IMMUNOASSAY SOLUTION



A robust market proven CLEIA analyzer utilising mono-test cartridges





FULLY AUTOMATED CHEMILUMINESCENT ENZYME IMMUNOASSAY SYSTEM

O1 CONSTANT THROUGHPUT

- 120 tests/hour independent of assay combinations or assay format order
- Each 30 seconds a result
- · Automated dilution



CARTRIDGE

- Interchangeable on all LUMIPULSE G series
- · Cost-transparent
- No wastage from on-board expiry
- · Simple operation



USER INTERFACE

- · Simple menu navigation
- · Limited hands-on time






PROVEN RELIABILITY
Meantime between failure
1125 days (figures April 2017)









NO SPECIMEN CARRY-OVER Instrument uses disposable sample tips equipped for level sensing








O 5 INCREASED WALK AWAY TIME
Load up to 504 tests and 100 samples




O 6 STABLE, SIMPLE CALIBRATORS

Most assays use 2-3 point calibrations and stable for 30 days



System architecture optimized for connection to automated sample tracking systems



Fujirebio Europe N.V.
Technologiepark 6 - 9052 Gent - Belgium T +32 (0)9 329 13 29 - F +32 (0)9 329 17 75 www.fujirebio-europe.com

