SIMPLIFY YOUR SINGLE CELL WORKFLOW LIKE NEVER BEFORE.

Meet Tecan's Uno Single Cell Dispenser™







SIMPLIFY YOUR SINGLE CELL WORKFLOW WITH TECAN'S UNO SINGLE CELL DISPENSER

Tecan's Uno Single Cell Dispenser is an automated benchtop instrument designed for ease of use, efficiency, and precision in single-cell workflows. Uno empowers a wide range of single-cell applications across the "-omics" landscape, including MS-based single-cell analysis, iPSC libraries, 3D cell research, cell-line development, and beyond.

With Uno, you can achieve accurate, viable, and precise single-cell dispensing, managing reagent volumes from picoliters to microliters within minutes.



Your needs at the core



Simple and accessible benchtop solution for labs of any size. Integrate intuitive software and compact hardware at a competitive price without compromising on performance.



Fast and precise dispensing saves reagents with picoliter-level accuracy. Isolate single cells in ~5 minutes¹ and dispense reagents in under three minutes across a 384-well plate.



Reliable dispensing process for reproducible results.

HP's proven microfluidic digital dispensing technology, a choice trusted by 1,000+ industry-leading customers.



Gentle single-cell dispensing to preserve cell integrity. Maintain a robust 90% cell viability¹ for downstream analysis. Works with a wide range of cells and reagents.

PRODUCT FEATURES

See what's possible for you

Uno Single Cell Dispenser Specification

Cell size range	10 μm to 25 μm
Cell types	Wide variety of cell types are tested
Qualified fluid classes	DMSO, Acetonitrile, aqueous solutions, aqueous 10 % trifluoroacetic acid (TFA), master mixes
Formats to dispense on	 12, 24, 48, 96, 384-well plates Customizable substrate, such as microscope slides Plate height 4-47 mm
Cell isolating technology	Microfluidic technology
Dispensing speed	 Cells: 384-well plate in ~5 minutes Reagents: 384-well plate in under three minutes
Single cell success rate	90% ²
Cell viability	90% 3
Software	Intuitive and easy to set-up protocolsAutomatically generates summary report, indicating wells with single cells
Overall instrument dimensions	47 cm x 38 cm x 23 cm
Instrument weight	4.8 kg (10.58 lb)

 $^{\rm 2)}$ Dispensing beads in 96-well plates $^{\rm 3)}$ At day zero



PRODUCT FEATURES

See what's possible for you

Uno Single Cell and Reagent Dispense Cassettes specification







C1a - Single Cell T1sf / T1 D1sf / D1 **Dispensehead Cassette Dispensehead Cassette Dispensehead Cassette Dispensing capability Cell Dispensing Reagent Dispensing Reagent Dispensing** Up to 400 single cells • T1sf: for surfactant-free • D1sf: for surfactant-free reagents reagents • T1: for surfactant-• D1: for surfactantrequiring reagents requiring reagents Minimum dispense volume 40 pl 11 pl 1 nl Fill cup size 20 µl 20 µl 200 µl

PRODUCT FEATURES

See what's possible for you

UnoControl Software

Designed with simplicity in mind.



Intuitive fluid and plate definition enables protocol set-up in minutes.

HeLa Cells #	Name HeLa Cells
CELL+S 200 cells/µL	Fluid group Cells 🗸
Lysis buffer	Fluid class Small (10-14 µm) 🗸
0.96 µL ©	Concentration 200 cells/µL 🗸
AqSF 10.0 mM	Dispense by # Count
Trypsin •	Source plate
0 µL © AqSF 10.0 mM	Source well
	Color 1 V
Create Diluted Fluids	Dispense 🌢 Yes

The protocol set-up guide includes a wide range of qualified fluids classes, and covers four distinct cell size ranges.



Guided steps ensure precise execution of experiments, minimizing handling errors.

General Dispense Summary Protocol (1998) Biological Statistical Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Statistical (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998) Control (1998)	Classes of Bulls and State St	
	2. SingleCellDispensing (3) xml	

Gain instant insights with auto-generated summary reports that detail dispensed reagent volumes, plates and cassettes used, single-cell isolation success rate, and cell concentration.

SIMPLIFY YOUR SINGLE CELL WORKFLOW LIKE NEVER BEFORE.



LEARN MORE



For research use only. Not for use in diagnostic procedures. Uno Single Cell Dispenser is a product of HP Inc. CA, USA.

Australia +61 3 9647 4100 Austria +43 62 46 89 330 Belgium +32 15 42 13 19 China +86 21 220 63 206 France +33 4 72 76 04 80 Germany +49 79 51 94 170 Italy +39 02 92 44 790 Japan +81 44 556 73 11 Netherlands +31 18 34 48 17 4 Nordic +46 8 750 39 40 Singapore +65 644 41 886 Spain +34 93 595 25 31 Switzerland +41 44 922 89 22 UK +44 118 9300 300 USA +1 919 361 5200 Other countries +41 44 922 81 11

Tecan Group Ltd. makes every effort to include accurate and up-to-date information within this publication, however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. All mentioned trademarks are protected by law. In general, the trademarks and designs referenced herein are trademarks, or registered trademarks, of Tecan Group Ltd., Männedorf, Switzerland. A complete list may be found at http://www.tecan.com/trademarks. Product names and company names that are not contained in the list but are noted herein may be the trademarks of their respective owners. For technical details and detailed procedures of the specifications provided in this document please contact your Tecan representative.

Tecan is in major countries a registered trademark of Tecan Group Ltd., Männedorf, Switzerland.

© 2023 Tecan Trading AG, Switzerland, all rights reserved.

www.tecan.com

